



ad'just®

LUMBAR CUSHION



Ad'just lumbar cushion

The solution to back pain in the car

A car seat should accommodate all body types. Whether you are short or tall, slim or stout, the seats are the same for everyone. They are not made to measure and do not provide care for your back. The right lumbar cushion adapts your seat to your body, keeps your body in line with the road and supports your lumbar arch.

THREE COMPONENTS

- **(A)** back piece
- **(B)** set of two ergonomic cushions (right and left)
- **(C)** inflatable lumbar support and pump

DIMENSIONS

See opposite

MATERIALS

- Quality materials selected for their technical capacity
- Fabric in contact with the user: "3D honeycomb polyester" mesh (dry / absorbent / aesthetic)
- Grip fabric: polyamide mesh (velvety, very sturdy appearance)
- Foam in the back piece: high resilience polyurethane (8 mm thick)
- Membrane / tube / pump: latex and PVC
- Ergonomic cushion padding: Expanded polystyrene microbeads
- Composition, washing and serial number labels

PACKAGING

Clear plastic zip bag.

Dimensions : 8 x 27 x 40 cm (approx..)

Weight: 600 g

MAINTENANCE

Hand washable. Ergonomic cushions and membrane with removable covers. Be careful to prevent contact between the velcro and the 3D fabric.

1- YEAR WARRANTY

Patent and trademark filed ®




contact@adjust.fr

www.adjust.fr

adjust sas : 92 route de Jallieu | 38080 Saint Marcel Bel Accueil | France

ad'just®

LUMBAR CUSHION



Ad'just lumbar cushion

The solution to back pain in the car

The ad'just seat cushion is recommended if your seat is tilted too far back.

The Ad'just cushion opens the angle between the torso and legs and releases tension in the lower back (detachable part at the tailbone to alleviate pressure). Stabilizing side cushions prevent the pelvis from tilting. Very harmful swaying left and right in turns !

Its firm, resilient foam absorbs vibrations.

It is provided with non-slip bands silicone. A back strap secures it to the lumbar cushion.

DIMENSIONS

See opposite

WEIGHT




A + B 580 gr

A 440 gr

B 140 gr

MATERIALS

- Quality materials selected for their technical capacity
- Fabric in contact with the user: "3D honeycomb polyester" mesh (dry / absorbent / aesthetic)
- Grip fabric: polyamide mesh (velvety, very sturdy appearance)
- Foam in the back piece: high resilience polyurethane (8 mm thick)
- Membrane / tube / pump: latex and PVC
- Ergonomic cushion padding: Expanded polystyrene microbeads
- Composition, washing and serial number labels



PACKAGING

Clear plastic bag.

Dimensions, h x l x L cm (approx) : 5x37x44

MAINTENANCE

Hand washable, removable covers.

1- YEAR WARRANTY



contact@adjust.fr

www.adjust.fr

adjust sas : 92 route de Jallieu | 38080 Saint Marcel Bel Accueil | France